


Cerium Bromide Detectors



Radiation Detectors



Introducing CeBr₃

Cerium Bromide is the latest development in room temperature, high resolution scintillation detectors and with a low intrinsic background, is a superior alternative to Lanthanum Halide based scintillators for applications including: homeland security, high resolution gamma spectrometry, medical imaging and geological exploration.

CeBr₃ has high light output, fast response, and shows very high energy and timing resolution. The detector's features include:

- Room temperature operation
- Same form factor as NaI(Tl) detectors
- Directly compatible with traditional scintillation detector electronics & multi-channel analysers
- Supplied in sealed metal housing

Above an energy of 200 keV the resolution is superior to NaI(Tl). It also shows good proportionality of response over the gamma energy range from 122 to 1274 keV*¹ and as cerium is an activator for the luminescence process, the crystal does not need any additional doping to aid scintillation, unlike Lanthanum Halides.

Furthermore CeBr₃ scintillation detectors do not suffer from the intrinsic background (self-activity) typical for Lanthanum Halide detectors like LaCl₃:Ce and LaBr₃:Ce.

In LaBr₃ for example, self-activity is primarily due to ¹³⁸La that emits conversion electrons and β-particles with energy of up to 1.7 MeV. The self-activity due to ¹³⁸La in LaBr₃ has an intrinsic count-rate of ~1.5 events/(cm³/sec). However, the self-activity of CeBr₃ (due to ¹⁴²Ce that emits β-particles with total energy of 4.5 MeV) is 4x10⁻⁴ events/(cm³/sec). Thus, the self-activity of CeBr₃ is about 3750 times lower than that in LaBr₃.^{*2}

While self-activity issues are not critical in PET (which relies on coincidence measurements), the negligible self-activity of CeBr₃ makes it more attractive in large radionuclide detection and imaging systems required for homeland security applications. In these systems, large detector volumes are required and the expected extrinsic count-rate can be extremely low which requires detectors with very low background.

CeBr₃ is well suited for applications requiring fast response, high count-rates, and good timing resolution. CeBr₃ scintillators can also be expected to provide accurate time-of-flight (TOF) information.

Uses

- Homeland Security
- Medical Imaging (PET)
- Gamma Spectrometry
- Single Photon Emission Computed Tomography (SPECT)
- Nuclear & Particle Physics Research
- X-Ray Diffraction (Very fast scintillation material)
- NDT - Non Destructive Testing
- Decommissioning & Environmental Monitoring
- Geological & Oil Exploration
- NORM Detection
- Food Contamination Monitoring

Properties


Density: 5.2 g / cc

Emission Wavelength: 380 nm


Decay Time: 17ns

Radiological Background: Negligible

Hygroscopic: Yes





4% Resolution at 662 keV



Low Intrinsic background

The Difference between NaI(Tl) and CeBr₃



	NaI(Tl)	CeBr ₃
Relative photoelectron yield	100	125
Energy (keV)	Typical energy resolution (FWHM) %	
30	18	22
60	12	15
81	11	13.5
122	9	11
356	8	5
662	6.5	4
1332	4.9	3
2600	4.0	2

References for this document:

*P.Schotanus - Scionix,

*K.S.Shah - New Scintillation Detectors for PET - 2010



• www.johncaunt.com • sales@johncaunt.com • +44 (0) 161 763 3334

© John Caunt Scientific Ltd 2017

